

NOW
AVAILABLE!

SOLARIS
FREEDOM TO LIVE

TheraFitTM

MANUFACTURED BY SOLARIS

ADJUSTABLE NIGHTWEAR GARMENTS

TheraFit's unique design gives medical professionals the ability to customize its compression level, simply by adjusting the laces. This provides patients with an effective and economical solution for nightly lymphoedema maintenance.

- Therapist-controlled gradient compression
- Customizable nightly therapy
- Simplifies home programs by eliminating multi-layered bandaging
- Comfortable and easy to care for
- Effective in maintaining the reductions achieved during MLD
- Designed by medical professionals




All Solaris products
are made in the USA

DISTRIBUTED BY:



Gooldscross, Cashel
Co. Tipperary
Ireland
Tel: 0504-42513
Fax: 0504-60098
office@tipptoes.com
www.tipptoes.com



TheraFit Upper Extremity Sizing Chart						
ITEM CODE	ITEM DESCRIPTION	ARM LENGTH	CIRCUMFERENCES (in Centimeters)			
			Palm	Wrist	Elbow	Axilla
TheraFit Arm						
TF-UE-AG-S-S-R/L	Short, Small, Right/Left	Up to 59.5	15 - 21	13.5 - 19	19 - 27.5	20 - 29
TF-UE-AG-S-M-R/L	Short, Medium, Right/Left	Up to 59.5	17 - 23	15.5 - 22	21.5 - 31	26 - 38
TF-UE-AG-S-L-R/L	Short, Large, Right/Left	Up to 59.5	19 - 27	17.5 - 25	24 - 35	32 - 47
TF-UE-AG-A-S-R/L	Avg., Small, Right/Left	Up to 68	15 - 21	13.5 - 19	19 - 27.5	20 - 29
TF-UE-AG-A-M-R/L	Avg., Medium, Right/Left	Up to 68	17 - 23	15.5 - 22	21.5 - 31	26 - 38
TF-UE-AG-A-L-R/L	Avg., Large, Right/Left	Up to 68	19 - 27	17.5 - 25	24 - 35	32 - 47
TF-UE-AG-T-S-R/L	Tall, Small, Right/Left	Up to 72	15 - 21	13.5 - 19	19 - 27.5	20 - 29
TF-UE-AG-T-M-R/L	Tall, Medium, Right/Left	Up to 72	17 - 23	15.5 - 22	21.5 - 31	26 - 38
TF-UE-AG-T-L-R/L	Tall, Large, Right/Left	Up to 72	19 - 27	17.5 - 25	24 - 35	32 - 47

To learn more or to order, please contact:



Gooldscross, Cashel
Co. Tipperary
Ireland
Tel: 0504-42513
Fax: 0504-60098
office@tipptoes.com
www.tipptoes.com