

Measuring Guidelines for Facial Units

Positioning Your Client When Measuring for Facial Garments

Tribute garments are most commonly worn while the client is sleeping, however, clients with facial swelling may be measured in a seated posture.

Measuring:

Please record the needed measurements in the designated areas on the Facial Order Form.

- A:** Measure the length from one ear tip to the other.
- B:** Measure the head circumference at the temple line or largest part.
- C:** Measure the distance from the outside crease of the right eye to the outside crease of the left eye.
- D:** Measure the length of the nose.
- E:** Measure the length at the midline from the nose to the upper lip.
- F:** Measure the length at the midline from the lower lip to the chin.
- G:** Measure the length from the sterno-clavicular notch to the chin. Follow the contour with chin in relaxed posture.
- H:** Measure from the right ear lobe to the right nostril.
- I:** Measure from the left ear lobe to the left nostril.
- J:** Measure the distance from the right TMJ to the chin.
- K:** Measure the distance from the left TMJ to the chin.
- L:** Measure the length from the right ear lobe to the HPS (High Point of Shoulder).
- M:** Measure the length from the left ear lobe to the HPS.
- N:** Measure the neck circumference.

Measuring Guidelines for Hand Units

Positioning Your Client When Measuring for Hand Garments

The Hand Order Forms are Right / Left specific.

Position your client's hand on a copy of the order form. Align the wrist crease to approximate the Center Wrist Line on the copy of the order form. Trace around the client's hand with a pen or pencil, keeping the tip of the pen/pencil close to the digits and hand. Have the client remove their hand from the copy of the Hand Order Form.

Measuring Circumferences:

Please record the needed measurements in the designated areas on the specific Hand Order Form.

- A:** Measure the wrist circumference over the styloid.
- B:** Measure the palm circumference across the metacarpal heads.
- C:** Measure the proximal circumference of the thumb.
- D:** Measure the proximal circumference of the index finger.
- E:** Measure the proximal circumference of the middle finger.
- F:** Measure the proximal circumference of the ring finger.
- G:** Measure the proximal circumference of the small finger.